

When ergonomic design
and intuition come together



ERAGON *axial*

Axial handle unites with modular system

ERAGON^{axial}


When ergonomic design and intuition come together

The **ERAGON modular** forceps and scissor system, already having a proven track record, has been expanded by the addition of a completely new axial handle. **ERAGON axial** combines all the benefits of the existing handle concept, such as wide finger contact surfaces, optimum operation of the rotation knob, precise grip, and secure hold ratchet lock. **ERAGON axial** offers fatigue-free working and precision with a new level of handling. See for yourself.



Ergonomic design gains upper hand

When the new handle was being created, simple handling and user-friendly operation were a top priority. The slightly angled contour permits a relaxed working position. The wide contact surfaces for the fingers facilitate fatigue-free working without the development of pressure points and therefore offer maximum surgeon comfort.



Turning with finesse

The optional knob attachment made of silicone makes operation even easier. A convenient accessory that is helpful for all hand sizes – small and large alike. The reusable attachment is simple to remove and easily reprocessed.

Safety under control

The intuitive, ratchet lock provides sensitive handling and a secure grip.

Skillfully combined

- **ERAGON axial** is fully compatible with all the components in the **ERAGON modular** Series. All sheaths and jaw inserts can be combined as required.
- The **"Click-it"** system provides the intuitive assembly and disassembly familiar for ERAGON.
- Make the most of the multi-faceted diversity of the successful **ERAGON modular** Series.

ERAGON^{axial}

When ergonomic design and intuition come together



ERAGON axial handle
with
ratchet lock 83930082



ERAGON axial handle
without lock, with
HF connector 83930083



ERAGON axial handle with
ratchet lock, with
HF connector 83930084



ERAGON axial handle
with ratchet lock,
not rotatable 83930085



Option:
Knob attachment
ERAGON axial,
autoclavable 8988

ERAGON^{modular}

All ERAGON modular components
are compatible with ERAGON axial.



WL = 450 mm

WL = 330 mm

WL = 260 mm

